

PROSCREEN™

Drugs of Abuse Test Cup




- ✚ **Made in the USA**
- ✚ **Zero Exposure to Donor Specimen - One-Step**
- ✚ **Test for up to 12 Drugs**
- ✚ **Donor Cannot Defeat Assay - Test Activated by Collector**
- ✚ **Built-in Specimen Validity Test Available**
- ✚ **Fewer False Positive Results from Common Medications**

Instructions for PROSCREEN™ Cup with On-Board Specimen Validity Test



Step 1

Have the donor provide a specimen in cup. (Collector retains test cup lid.)



Step 2

Verify that temperature is within range on the temperature strip (90° - 100° F).




Step 3

Open foil pouch containing test lid. Have the donor screw the cap onto cup until tight.



Step 4

Tilt Cup on side to rest on tabs. Test will run automatically.



Step 5

On cups equipped with specimen validity test, read specimen validity results at the times specified on the color chart provided. If test pad colors are outside the normal range, open a new cup, collect a second specimen, and retest.



Step 6

Read the drug test results when the control lines are clearly visible.

Interpreting the Results



NEGATIVE

A *Negative* result is indicated by the presence of a test or "T" line for each designated drug.



NEGATIVE

The presence of a very light "T" line indicates a *Negative* result.



POSITIVE

A *Positive* result is indicated by the presence of a control "C" line and the absence of a "T" line. Wait 5 minutes to read positive results.



INVALID

An *Invalid* result is indicated when the "C" line is completely missing from one or more test windows. If this happens, run another test.